



mustangTM

The Mustang is a walking aid for disabled children and youngsters, who need supplementary assistance with their daily walking



The Mustang provides excellent and snug support for the upper body whilst also keeping the child's pelvis gently supported. The upright, slightly forward leaning position stimulates the walking function and encourages the child to participate in daily activities.

The Mustang can be used by e.g. CP children and others with lack of balance and impairment of their motor function. The Mustang is developed as a basic anterior walker but can also be used as a posterior gait trainer, depending on the child's level of ability.

The Mustang is available in four sizes and suitable for children and young people.



14-60 cm




37-125 cm



max.130 cm



max.80 kg



Motivating the child to walk

The user is placed at the centre of the Mustang which provides core stability and enables the user to focus energy on walking



Develop a more natural walking pattern

The Mustang has easy and steady angle adjustable centre bar. A slight forward leaning position, combined with built-in suspension, stimulates the natural walking pattern



Easy to access

The Mustang has been constructed with accessibility in mind. The U-shaped frame and wrap around support, allows for easy access and transfer



Easy to store and carry

The foldable frame makes storage and transport of the Mustang easy



Standard

As standard the Mustang is delivered with lockable wheels. The Mustang is available in 4 sizes



Size 1 is available in yellow and with 4" or 8" wheels



Size 2 is available in red, white or erbium and with 4" or 8" wheels

Standard




Size 3 is available in red, white or erbium and with 8" wheels



Size 4 is available in erbium and with 8" wheels

Frame colour

Available in 4 colours



Yellow = Size 1



Red = size 2+3



White = size 2+3



Erbium = size 2+3+4

Wheels

The Mustang is delivered as standard with various adjustment possibilities on all four wheels



Directional stabilizer



Tension adjustable wheels



Wheel brake on all 4 wheels



Lock




Release

Wheel stops

The unique on/off anti-reverse wheel stops are very easy to lock and unlock



nothing compares to a smile



Measurements

	Size 1	Size 2	Size 3	Size 4
Width (G)	52 cm (20 ³ / ₄ "	58 cm (22 ¹ / ₂ "	65 cm (25 ¹ / ₂ "	70 cm (27 ¹ / ₂ "
Width between wheels (G1)	38 cm (14 ³ / ₄ "	45 cm (17 ¹ / ₂ "	50 cm (19 ¹ / ₂ "	56 cm (22"
Length with 4" wheels (H)	75 cm (29 ¹ / ₄ "	83 cm (32 ¹ / ₂ "	-	-
Length with 8" wheels (H)	82 cm (32"	90 cm (35"	100 cm (39"	110 cm (43"
Height with 4" wheels (I)	40-65 cm (14 ¹ / ₄ -25 ¹ / ₄ "	60-85 cm (23 ¹ / ₂ -33 ¹ / ₂ "	-	-
Height with 8" wheels (I)	49-65 cm (19-25 ¹ / ₄ "	60-85 cm (23 ¹ / ₂ -33 ¹ / ₂ "	70-105 cm (27 ¹ / ₂ -41"	87-130 cm (34-50 ³ / ₄ "
Seat height above floor, with 4" wheels (seat = accessory)	17-35 cm (6 ¹ / ₂ -13 ¹ / ₂ "	24-41 cm (9 ¹ / ₂ -16"	-	-
Seat height above floor, with 8" wheels (seat = accessory)	26 - 44 cm (10-17 ¹ / ₄ "	33-50 cm (12 ³ / ₄ -19 ¹ / ₂ "	45-68 cm (17 ¹ / ₂ -26 ¹ / ₂ "	62-84 cm (24-32 ³ / ₄ "
Width between handles (handles = accessory)	14-35 cm (5 ¹ / ₂ -13 ¹ / ₂ "	14-35 cm (5 ¹ / ₂ -13 ¹ / ₂ "	14-60 cm (5 ¹ / ₂ -23 ¹ / ₂ "	14-60 cm (5 ¹ / ₂ -23 ¹ / ₂ "
Angle of the centre bar (a1)	30°	30°	0-25°	0-25°
Turning radius	46 cm (18"	55 cm (21 ¹ / ₂ "	60 cm (23 ¹ / ₂ "	65 cm (25 ¹ / ₄ "
User height, floor to axil	40-65 cm (14 ¹ / ₄ -25 ¹ / ₄ "	60-85 cm (23 ¹ / ₂ -33 ¹ / ₂ "	70-105 cm (27 ¹ / ₂ -41"	87-130 cm (34-50 ³ / ₄ "
Weight with 4" wheels, standard without accessories	10 kg (22 lb)	11 kg (24 lb)	-	-
Weight with 8" wheels, standard without accessories	13 kg (28,5 lb)	14 kg (31 lb)	18 kg (40 lb)	19 kg (42 lb)
Max. user weight / load	30 kg (66 lb)	40 kg (88 lb)	60 kg (132 lb)	80 kg (176 lb)

ACCESSORIES



Chest support
8698x0
For support of the upper body

size	measurements
S	O: 47-64 cm
M	O: 67-78 cm
L	O: 87-114 cm
O:	inner circumference



Cover
8698x0-x60
For added comfort

size	measurements
1	H:8½ x W:25 cm
2	H:12 x W:38 cm
3	H:15 x W:49 cm



Extension piece
869225
To eliminate the gap when width adjusting the chest support

size	measurements
one	H: 11 x W: 5½ cm



Hand bar
86936x
Provide hand control For anterior use

size	measurements
S	H:11 x W:21 cm
L	H:11 x W:32 cm



Abdominal support
869260
Provides additional anterior pelvic control

size	measurements
one	H:13 x W:19 cm



Arm rests, anterior
86935x
Provide forearm support. Helps the child to stabilise the shoulder girdle

Two sizes



Arm prompts, anterior or posterior
86933x
Forearm support to help stabilise the shoulder girdle

Three sizes



Padded forearm positioner
869370
For the better grip and weight bearing

One size



Sacrum support
86925x
Helps to stabilise the pelvis. Requires seat 8691X0

size	measurements
1	H:18 cm (from seat) D:13-21 cm (from seat)
2	H:20½ cm (from seat) D:16½-24½ cm (f.seat)



Hip support, small
869220 or 869221
Helps to stabilise the pelvis. Requires seat 8691X0

size	measurements
one	W:18½-25½ cm



Hip support, large
869230 or 869231
Helps to stabilise the pelvis. Requires seat 8691X0

size	measurements
one	W:25½-32½ cm



Connecting rod for head support system
869226
Mounted between hip and head support

One size



ACCESSORIES



Seat
8691x0
Provides a resting position that a child can use when fatigued. Used for mounting sacrum or hip supports

size	measurements
1	H:17½ x D:25½ cm
2	H:17½ x D:30½ cm



Sling seat
86957x
Provides partial weight bearing support

size	measurements
1	W:18 x L:24½ cm
2	W:18 x L:33 cm



Leg pad
869470
Provides knee protection

One size



Guide bar
869510
Allows caregiver to direct and assist the child's movements.
For Mustang size 1 and 2

size	measurements
one	L: 70 cm



Storage pouch
869560
For transporting toys etc. 0.5 kg weight limit.
Only for size 1+2

size	measurements
one	H:14½ x W:13 cm



Storage bag
869565
Only for size 3+4. 3 kg weight limit

size	measurements
one	H: 30 x W: 39½ cm



Leg separator
8694x0
Helps to prevent feet crossing midline and so maintain a neutral walking pattern.

Four sizes



Ankle prompts
8694xx
Guidance for correct walking pattern

size	measurements
1	L: 52½ cm
2	L: 60 cm
3	L: 70½ cm
4	L: 80½ cm

